





















## MINI FRAGMENTS IMPLANTS (SS 316L)

### 1.5mm Screws & Plates









<p><b>T001</b> Code No. T001.04 T001.05 T001.06 T001.07 T001.08 T001.10</p>	<p><b>1.5mm Micro Cortex Screws</b> <u>Total Length in mm</u> 4mm 5mm 6mm 7mm 8mm 10mm</p>	<p>Self Taping Threaded Diameter 1.5mm Drill Bit for Threading Hole 1.1mm Core Diameter 1.1mm Head Diameter 2.5mm Cruciform Recess Head</p>	
<p><b>T002</b> Code No. T002.02 T002.03 T002.04 T005.05</p>	<p><b>1.5mm Micro Plate with Center Space</b> <u>Total Hole</u> 2hole 3hole 4hole 5hole</p>		<p>Plate Thickness 0.8mm For 1.5mm Screws With Middle Hole spacing</p> 
<p><b>T003</b> Code No. T003.04 T003.05 T003.06 T003.20</p>	<p><b>1.5mm Micro Plate (Tit)</b> <u>Total Hole</u> 4hole 5hole 6hole 20hole</p>	<p>Plate Thickness 0.8mm For 1.5mm Screws</p>	
<p><b>T004</b> Code No. T004.04 T004.06 T004.08</p>	<p><b>1.5mm Micro Orbital Plate with center Space(Tit)</b> <u>Total Hole</u> 4hole 6hole 8hole</p>		<p>Plate Thickness 0.8mm For 1.5mm Screws With Middle Hole spacing</p> 
<p><b>T005</b> Code No. T005.08 T005.10 T005.12</p>	<p><b>1.5mm Micro Orbital Plate (Tit)</b> <u>Total Hole</u> 8hole 10hole 12hole</p>	<p>Plate Thickness 0.8mm For 1.5mm Screws</p>	
<p><b>T006</b> Code No. T006.04L T006.04R</p>	<p><b>1.5mm Micro L Plate Angle 2hole (Tit)</b> <u>Total Hole</u> 2hole x 2hole (Left) 2hole x 2hole (Right)</p>		<p>Plate Thickness 0.8mm For 1.5mm Screws With Middle Hole spacing</p> 
<p><b>T007</b> Code No. T007.06L T007.06R</p>	<p><b>1.5mm Micro L Plate Angle 3hole (Tit)</b> <u>Total Hole</u> 3hole x 3hole (Left) 3hole x 3hole (Right)</p>	<p>Plate Thickness 0.8mm For 1.5mm Screws With Middle Hole spacing</p>	
<p><b>T008</b> Code No. T008.04</p>	<p><u>Total Hole</u> 2hole x 2hole</p>		<p>Plate Thickness 0.8mm For 1.5mm Screws With Middle Hole spacing</p> 
<p><b>T009</b> Code No. T009.05</p>	<p><b>1.5mm Micro X Plate (Tit)</b> <u>Total Hole</u> 5hole</p>	<p>Plate Thickness 0.8mm For 1.5mm Screws With Middle Hole spacing</p>	
<p><b>T010</b> Code No. T010.05</p>	<p><b>1.5mm Micro Y Plate (Tit)</b> <u>Total Hole</u> 5hole</p>		<p>Plate Thickness 0.8mm For 1.5mm Screws With Middle Hole spacing</p> 





**2.0mm Screws & Plates**

<p><b>T025</b> <u>Code No.</u> T025.05 T025.06 T025.07 T025.08 T025.09 T025.10 T025.12</p>	<p><b>2mm Cortex Screws (Tit)</b> <u>Total Length in mm</u> 5mm 6mm 7mm 8mm 9mm 10mm 12mm</p>	<p>Self Taping Threaded Diameter 2.0mm Drill Bit for Threading Hole 1.5mm Core Diameter 1.4mm Head Diameter 3.0mm Cruciform Recess Head</p>	
<p><b>T026</b> <u>Code No.</u> T026.02 T026.03 T026.04 T026.05 T026.06</p>	<p><b>2mm Mandible Plate with Center Space (Tit)</b> <u>Total Hole</u> 2hole 16mm 3hole 22mm 4hole 28mm 5hole 34mm 6hole 40mm</p>	<p>Thickness 1.0mm Width 5.0mm Hole Spacing 6.0mm Middle Hole spacing 11.0mm For 2.0mm Screws</p>	
<p><b>T027</b> <u>Code No.</u> T027.04 T027.05 T027.06 T027.08 T027.10 T027.20</p>	<p><b>2mm Mandible Plate (Tit)</b> <u>Total Hole</u> 4hole 23mm 5hole 29mm 6hole 35mm 8hole 47mm 10hole 59mm 20hole 119mm</p>	<p>Thickness 1.0mm Width 5.0mm Hole Spacing 6.0mm For 2.0mm Screws</p>	
<p><b>T028</b> <u>Code No.</u> T028.04 T028.06 T028.08</p>	<p><b>2mm Orbital Plate with Center Space(Tit)</b> <u>Total Hole</u> 4hole 6hole 8hole</p>	<p>Thickness 1.0mm Width 5.0mm Hole Spacing 6.0mm Middle Hole spacing 11.0mm For 2.0mm Screws</p>	
<p><b>T029</b> <u>Code No.</u> T029.08 T029.10 T029.12</p>	<p><b>2mm Orbital Plate (Tit)</b> <u>Total Hole</u> 8hole 10hole 12hole</p>	<p>Thickness 1.0mm Width 5.0mm Hole Spacing 6.0mm For 2.0mm Screws</p>	
<p><b>T030</b> <u>Code No.</u> T030.04L T030.04R T030.05L T030.05R T030.06L T030.06R</p>	<p><b>2mm L Plate Angle 2hole (Tit)</b> <u>Total Hole</u> 2hole x 2hole (Left) 2hole x 2hole (Right) 2hole x 3hole (Left) 2hole x 3hole (Right) 2hole x 4hole (Left) 2hole x 4hole (Right)</p>	<p>Thickness 1.0mm Width 5.0mm Hole Spacing 6.0mm Middle Hole spacing 11.0mm For 2.0mm Screws</p>	
<p><b>T031</b> <u>Code No.</u> T031.06L T031.06R</p>	<p><b>2mm L Plate Angle 3hole</b> <u>Total Hole</u> 3hole x 3hole (Left) 3hole x 3hole (Right)</p>	<p>Thickness 1.0mm Width 5.0mm Hole Spacing 6.0mm Middle Hole spacing 11.0mm For 2.0mm Screws</p>	
<p><b>T032</b> <u>Code No.</u> T032.04 T032.05 T032.06</p>	<p><b>2mm T Plate Head 2 Hole</b> <u>Total Hole</u> 2hole x 2hole 2hole x 3hole 2hole x 4hole</p>	<p>Thickness 1.0mm Width 5.0mm Hole Spacing 6.0mm Middle Hole spacing 11.0mm For 2.0mm Screws</p>	
<p><b>T036</b> <u>Code No.</u> T036.05</p>	<p><b>2mm X Plate</b> <u>Total Hole</u> 5hole</p>	<p>Thickness 1.0mm Width 5.0mm Middle Hole spacing 11.0mm For 2.0mm Screws</p>	
<p><b>T037</b> <u>Code No.</u> T037.04 T037.05</p>	<p><b>2mm Y Plate</b> <u>Total Hole</u> 4hole 5hole</p>	<p>Thickness 1.0mm Width 5.0mm Hole Spacing 6.0mm Middle Hole spacing 11.0mm For 2.0mm Screws</p>	




<b>T039</b> Code No. T039.04L T039.04R	<b>2mm Z Plate</b>  Total Hole 4hole (Left) 4hole (Right)	Thickness 1.0mm Width 5.0mm Hole Spacing 6.0mm Middle Hole spacing 11.0mm For 2.0mm Screws	 Left Right
---	---	--	---

### 2.4mm Screws & Plates

<b>T051</b> Code No. T051.06 T051.07 T051.08 T051.09 T051.10 T051.12 T051.14	<b>2.4mm Cortex Screws</b>  Total Length in mm 6mm 7mm 8mm 9mm 10mm 12mm 14mm	Self Taping Threaded Diameter 2.4mm Drill Bit for Threading Hole 1.8mm Core Diameter 1.7mm Head Diameter 4.0mm Cruciform Recess Head	
<b>T052</b> Code No. T052.02 T052.03 T052.04 T052.05 T052.06	<b>2.4mm Straight Plate with Center Space</b>  Total Hole 2hole 19mm 3hole 25mm 4hole 31mm 5hole 37mm 6hole 43mm	Thickness 1.0mm Width 6.0mm Hole Spacing 6.0mm Middle Hole spacing 13.5mm For 2.4mm Screws	
<b>T053</b> Code No. T053.03 T053.04 T056.05 T053.06 T053.20	<b>2.4mm Straight Plate</b>  Total Hole 3hole 18mm 4hole 24mm 5hole 30mm 6hole 36mm 20hole 120mm	Thickness 1.0mm Width 6.0mm Hole Spacing 6.0mm For 2.4mm Screws	
<b>T054</b> Code No. T054.06 T054.08	<b>2.4mm Orbital Plate With Center Space</b>  Total Hole 6hole 8hole	Thickness 1.0mm Width 6.0mm Hole Spacing 6.0mm Middle Hole spacing 13.5mm For 2.4mm Screws	
<b>T055</b> Code No. T055.04L T055.04R	<b>2.4mm L Plate Angle 2hole</b>  Total Hole 2hole x 2hole (Left) 2hole x 2hole (Right)	Thickness 1.0mm Width 6.0mm Hole Spacing 6.0mm Middle Hole spacing 13.5mm For 2.4mm Screws	 Left Right
<b>T056</b> Code No. T056.04 T056.05 T056.06	<b>2.4mm T Plate Head 2hole</b>  Total Hole 2hole x 2hole 2hole x 3hole 2hole x 4hole	Thickness 1.0mm Width 6.0mm Hole Spacing 6.0mm Middle Hole spacing 13.5mm For 2.4mm Screws	
<b>T057</b> Code No. T057.05	<b>2.4mm X Plate</b>  Total Hole 5hole	Thickness 1.0mm Width 6.0mm Middle Hole spacing 13.5mm For 2.4mm Screws	
<b>T058</b> Code No. T058.04 T058.05	<b>2.4mm Y Plate</b>  Total Hole 4hole 5hole	Thickness 1.0mm Width 6.0mm Hole Spacing 6.0mm Middle Hole spacing 13.5mm For 2.4mm Screws	

<b>T059</b>		<b>2.4mm Double Y Plate</b>		Thickness	1.0mm		
<u>Code No.</u>	<u>Total Hole</u>	Width	6.0mm	Hole Spacing	6.0mm		
T059.08	8hole	Middle Hole spacing	13.5mm	For 2.4mm Screws			
<b>T060</b>		<b>2.4mm Z Plate</b>		Thickness	1.0mm		
<u>Code No.</u>	<u>Total Hole</u>	Width	6.0mm	Hole Spacing	6.0mm		
T060.04L	4hole (Left)	Middle Hole spacing	13.5mm	For 2.4mm Screws			
T060.04R	4hole (Right)						
<b>T061</b>		<b>Mini Straight Reconstruction Plate 2.4mm</b>		Thickness	1.0mm		
<u>Code No.</u>	<u>Total Hole</u>	Width	6.0mm	Hole Spacing	6.0mm		
T061.04	4hole 24mm	For 2.4mm Screws					
T061.06	6hole 36mm						
T061.08	8hole 48mm						
T061.20	20hole 120mm						
<b>T062</b>		<b>Angled Reconstruction Plate 2.4mm</b>		Thickness	1.0mm		
<u>Code No.</u>	<u>Total Hole</u>	Width	6.0mm	Hole Spacing	6.0mm		
T062.03	3+3hole	For 2.4mm Screw					
T062.04	4+4hole						
T062.16L	4+16hole(Left)						
T062.16R	4+16hole(Right)						

## 2.7mm Screws & Plates

<b>076</b>		<b>T076</b>		<b>2.7mm Cortex Screws</b>				
<u>Code No.</u>	<u>Code No.</u>	<u>Total Length in mm</u>	Thickness	2.7mm	Drill Bit for Threading Hole		2.0mm	
076.08	T076.08	8mm	Width	8.0mm	Drill Bit for Gliding Hole		2.7mm	
076.10	T076.10	10mm	Hole Spacing	8.0mm	Core Diameter		1.9mm	
076.12	T076.12	12mm	Head Diameter				5.0mm	
076.14	T076.14	14mm	Hexagonal Head Socket				2.5mm	
076.16	T076.16	16mm						
076.18	T076.18	18mm						
076.20	T076.20	20mm						
076.22	T076.22	22mm						
076.24	T076.24	24mm						
076.26	T076.26	26mm						
076.28	T076.28	28mm						
076.30	T076.30	30mm						
<b>077</b>		<b>T077</b>		<b>Mini Straight Reconstruction Plate 2.7mm</b>				
<u>Code No.</u>	<u>Code No.</u>	<u>Total Length in mm</u>	Thickness	2.0mm				
077.04	T077.04	4hole 32mm	Width	8.0mm				
077.05	T077.05	5hole 40 mm	Hole Spacing	8.0mm				
077.06	T077.06	6hole 48mm	For Use 2.7mm Screws					
077.07	T077.07	7hole 56mm						
077.08	T077.08	8hole 64mm						
077.09	T077.09	9hole 72mm						
077.10	T077.10	10hole 80mm						
<b>078</b>		<b>T078</b>		<b>Mini DCP 2.7mm</b>				
<u>Code No.</u>	<u>Code No.</u>	<u>Total Length in mm</u>	Thickness	2.0mm				
078.02	T078.02	2hole 20mm	Width	8.0mm				
078.03	T078.03	3hole 28mm	Hole Spacing	8.0mm				
078.04	T078.04	4hole 36mm	For Use 2.7mm Screws					
078.05	T078.05	5hole 44 mm						
078.06	T078.06	6hole 52mm						
078.07	T078.07	7hole 60mm						
078.08	T078.08	8hole 68mm						

This document was created with Win2PDF available at <http://www.daneprairie.com>.  
The unregistered version of Win2PDF is for evaluation or non-commercial use only.