

## CANNULATED SCREW & INSTRUMENTS

### 4.0mm Cannulated Screws Instruments Set

Guide Wire 1.25mm	10 No.
Drill Bit Cannulated 2.7mm	1 No.
Parallel Drill Guide 1.25mm	1 No.
Direct Measuring Device for 1.25mm Guide Wire	1 No.
Tap Cannulated for 4.0mm Screws	1 No.
Screw Driver Cannulated 2.5mm Tip	1 No.
Mini Screw Forceps Self-Holding	1 No.

### 4.0mm Cannulated Screws Implant Set

4.0mm Cannulated Screw (Short Thread) 10mm, 12mm, 14mm, 16mm, 18mm, 20mm, 22mm, 24mm, 26mm, 28mm, 30mm, 32mm, 34mm, 36mm, 38mm, 40mm, 42mm, 44mm, 46mm, 48mm, 50mm,	Each 3 No.
4.0mm Cannulated Screw (Fully Thread) 10mm, 12mm, 14mm, 16mm, 18mm, 20mm, 22mm, 24mm, 26mm, 28mm, 30mm, 32mm, 34mm, 36mm, 38mm, 40mm, 42mm, 44mm, 46mm, 48mm, 50mm,	Each 3 No.
Washer 7.0/3.6mm Dia for small Screws.	6 No.

### 4.0mm Cannulated Screw

306 Code No.	T306 Code No.	4.0mm Cannulated Screw (Short Thread) Total Length in mm / Thread Length		
306.10	T306.10	10mm	4mm	
306.12	T306.12	12mm	4mm	
306.14	T306.14	14mm	5mm	
306.16	T306.16	16mm	5mm	
306.18	T306.18	18mm	6mm	
306.20	T306.20	20mm	7mm	
306.22	T306.22	22mm	7mm	Thread Diameter 4.0mm
306.24	T306.24	24mm	8mm	Drill Bit for Threading Hole 2.7mm
306.26	T306.26	26mm	8mm	Shaft Diameter 2.7mm
306.28	T306.28	28mm	9mm	Core Diameter 2.5mm
306.30	T306.30	30mm	10mm	Hexagonal Socket Width 2.5mm
306.32	T306.32	32mm	11mm	Head Diameter 6.0mm
306.34	T306.34	34mm	11mm	Cannulation 1.35mm
306.36	T306.36	36mm	12mm	For Use Guide wire 1.2mm
306.38	T306.38	38mm	12mm	
306.40	T306.40	40mm	13mm	
306.45	T306.45	45mm	14mm	
306.50	T306.50	50mm	16mm	
306.55	T306.55	55mm	18mm	
306.60	T306.60	60mm	20mm	
306.65	T306.65	65mm	22mm	
306.70	T306.70	70mm	24mm	



### 6.5mm Cannulated Screws Instruments Set

Guide Wire 2.0mm	10 No.
Drill Bit Cannulated 4.5mm	1 No.
Parallel Drill Guide 2.0mm	1 No.
Direct Measuring Device for 2.0mm Guide Wire	1 No.
Tap Cannulated for 6.5mm Screws	1 No.
Screw Driver Cannulated For 4.5mm Screws	1 No.
Screw Forceps Self-Holding	1 No.

### 6.5mm Cannulated Screws Implant Set

6.5mm Cannulated Screw (16mm Thread) 25mm, 30mm, 35mm, 40mm, 45mm, 50mm, 55mm, 60mm, 65mm, 70mm, 75mm	Each 2 No.
80mm, 85mm, 90mm,	Each 4 No.
95mm, 100mm, 105mm., 110mm, 115mm, 120mm	Each 2 No.
6.5mm Cannulated Screw (32mm Thread) 45mm, 50mm, 55mm, 60mm, 65mm, 70mm, 75mm, 80mm, 85mm, 90mm,	Each 2 No.
95mm, 100mm, 105mm., 110mm, 115mm, 120mm	Each 1 No.
Washer 13.0/6.6mm for Large Screws	6 No.

## 6.5mm Cannulated Screws

<b>310 T310</b>		<b>6.5mm Cannulated Screw (16mm Thread)</b>			
<u>Code No.</u>	<u>Code No.</u>	<u>Total Length in mm</u>			
310.25	T310.25	25mm			
310.30	T310.30	30mm			
310.35	T310.35	35mm			
310.40	T310.40	40mm			
310.45	T310.45	45mm			
310.50	T310.50	50mm			
310.55	T310.55	55mm			
310.60	T310.60	60mm			
310.65	T310.65	65mm			
310.70	T310.70	70mm			
310.75	T310.75	75mm			
310.80	T310.80	80mm			
310.85	T310.85	85mm			
310.90	T310.90	90mm			
310.95	T310.95	95mm			
310.100	T310.100	100mm			
310.105	T310.105	105mm			
310.110	T310.110	110mm			
310.115	T310.115	115mm			
310.120	T310.120	120mm			

<b>311 T311</b>		<b>6.5mm Cannulated Screw (32mm Thread)</b>			
<u>Code No.</u>	<u>Code No.</u>	<u>Total Length in mm</u>			
311.45	T311.45	45mm			
311.50	T311.50	50mm			
311.55	T311.55	55mm			
311.60	T311.60	60mm			
311.65	T311.65	65mm			
311.70	T311.70	70mm			
311.75	T311.75	75mm			
311.80	T311.80	80mm			
311.85	T311.85	85mm			
311.90	T311.90	90mm			
311.95	T311.95	95mm			
311.10	T311.10	100mm			
311.105	T311.105	105mm			
311.110	T311.110	110mm			
311.115	T311.115	115mm			
311.120	T311.120	120mm			

Thread Diameter	6.5mm
Thread Length	16.0mm
Drill Bit for Threading Hole	4.5mm
Shaft Diameter	4.5mm
Core Diameter	4.5mm
Hexagonal Socket Width	3.5mm
Head Diameter	8.0mm
Cannulation	2.1mm
For Use Guide wire 1.8mm	

Thread Diameter	6.5mm
Thread Length	32.0mm
Drill Bit for Threading Hole	4.5mm
Shaft Diameter	4.5mm
Core Diameter	4.5mm
Hexagonal Socket Width	3.5mm
Head Diameter	8.0mm
Cannulation	2.1mm
For Use guide wire 1.8mm	

## 7.0mm Cannulated Screws Instruments Set


Guide Wire 2.0mm	10 No.
Drill Bit Cannulated 4.5mm	1 No.
Parallel Drill Guide 2.0mm	1 No.
Direct Measuring Device for 2.0mm Guide Wire	1 No.
Tap Cannulated for 7.0mm Screws	1 No.
Screw Driver Cannulated For 4.5mm Screws	1 No.
Screw Forceps Self-Holding	1 No.

## 7.0mm Cannulated Screws Implant Set

7.0mm Cannulated Screw (16mm Thread) 30mm, 35mm, 40mm, 45mm, 50mm, 55mm, 60mm, 65mm, 70mm, 75mm	Each 2 No.
80mm, 85mm, 90mm,	Each 4 No.
95mm, 100mm, 105mm., 110mm, 115mm, 120mm	Each 2 No.
7.0mm Cannulated Screw (32mm Thread) 45mm, 50mm, 55mm, 60mm, 65mm, 70mm, 75mm, 80mm, 85mm, 90mm,	Each 2 No.
95mm, 100mm, 105mm., 110mm, 115mm, 120mm	Each 1 No.
Washer 13.0/6.6mm for Large Screws	6 No.

## 7.0mm Cannulated Screws

<b>313</b>		<b>T313</b>	<b>7.0mm Cannulated Screw (16mm Thread)</b>	
<u>Code No.</u>	<u>Code No.</u>	<u>Code No.</u>	<u>Total Length in mm</u>	
313.30	T313.30		30mm	
313.35	T313.35		35mm	
313.40	T313.40		40mm	
313.45	T313.45		45mm	
313.50	T313.50		50mm	
313.55	T313.55		55mm	
313.60	T313.60		60mm	
313.65	T313.65		65mm	
313.70	T313.70		70mm	
313.75	T313.75		75mm	
313.80	T313.80		80mm	
313.85	T313.85		85mm	
313.90	T313.90		90mm	
313.95	T313.95		95mm	
313.100	T313.100		100mm	
313.105	T313.105		105mm	
313.110	T313.110		110mm	
313.115	T313.115		115mm	
313.120	T313.120		120mm	



Thread Diameter	7.0mm
Thread Length	16.0mm
Drill Bit for Threading Hole	4.5mm
Shaft Diameter	4.5mm
Core Diameter	4.5mm
Hexagonal Socket Width	3.5mm
Head Diameter	8.0mm
Cannulation	2.1mm
For Use Guide wire	1.8mm


<b>314</b>		<b>T314</b>	<b>6.5mm Cannulated Screw (32mm Thread)</b>	
<u>Code No.</u>	<u>Code No.</u>	<u>Code No.</u>	<u>Total Length in mm</u>	
314.40	T314.40		45mm	
314.45	T314.45		50mm	
314.50	T314.50		55mm	
314.55	T314.55		60mm	
314.60	T314.60		65mm	
314.65	T314.65		70mm	
314.70	T314.70		75mm	
314.75	T314.75		80mm	
314.80	T314.80		85mm	
314.85	T314.85		90mm	
314.90	T314.90		95mm	
314.95	T314.95		100mm	
314.100	T314.100		105mm	
314.105	T314.105		110mm	
314.110	T314.110		115mm	
314.115	T314.115		120mm	
314.120	T314.120			

Thread Diameter	6.5mm
Thread Length	32.0mm
Drill Bit for Threading Hole	4.5mm
Shaft Diameter	4.5mm
Core Diameter	4.5mm
Hexagonal Socket Width	3.5mm
Head Diameter	8.0mm
Cannulation	2.1mm
For Use guide wire	1.8mm




## Cannulated Interference Screw

<b>T315</b>	<b>Cannulated Interference Screw (ACL Screw) Titanium</b>		
<u>Length in mm</u>	<b>7mm Code No. T315A</b>	<b>8mm Code No. T315B</b>	<b>9mm Code No. T315C</b>
20mm	T315A.20	T315B.20	T315C.20
25mm	T315A.25	T315B.25	T315C.25
30mm	T315A.30	T315B.30	T315C.30
35mm	T315A.35	T315B.35	T315C.35



## Herbert Screw 3.5mm Cannulated

<b>T316</b>	<b>Herbert Screw 3.5mm Cannulated Titanium</b>	
<u>Code No.</u>	<u>Total Length in mm</u>	
T316.14	14mm	
T316.16	16mm	
T316.18	18mm	
T316.20	20mm	
T316.22	22mm	
T316.24	24mm	
T316.26	26mm	
T316.28	28mm	
T316.30	30mm	



This document was created with Win2PDF available at <http://www.daneprairie.com>.  
The unregistered version of Win2PDF is for evaluation or non-commercial use only.